



COVID-19/Flu A&B/RSV Ag Combo Rapid Test Cassette (Swab)

► Product Information

The COVID-19/Flu A&B/RSV Ag Combo Rapid Test Cassette (Swab) is an in vitro immunochromatographic assay for the qualitative and differential detection of nucleocapsid protein antigen from SARS-CoV-2, influenza A, influenza B and Respiratory Syncytial Virus (RSV) in nasal swab specimens directly from individuals within the first ten days of symptom onset. It is intended to aid in the rapid diagnosis of SARS-CoV-2, influenza A, RSV infections.

Benefits

- Rapid testing for SARS-CoV-2, Influenza A&B, RSV antigen within 15 minutes
- Facilitates patient treatment decisions quickly
- Simple, time-saving procedure
- All necessary reagents provided & no equipment needed
- High sensitivity and specificity

Specification


Information	Detail
Time to result	15 minutes
Storage	2-30°C
Shelf life	24 months
Specimen type	Nasal Swab

Contents

- 20 Test cassettes
- 20 Sterile swabs
- 20 Extraction tubes with buffer and tips
- 2 Workstations
- 1 Package insert

Test Procedure & Interpretation

Specimen Preparation



Step 1

Using the sterile swab provided in the kit, carefully insert the swab into one nostril of the patient. The swab tip should be inserted up to 2-4 cm until resistance is met.



Step 2

Roll the swab 5 times along the mucosa inside the nostril to ensure that both mucus and cells are collected.



Step 3


Using the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities.



Step 4

Withdraw the swab from the nasal cavity. The specimen is now ready for preparation using the extraction buffer provided in the test kit.

Sample Preparation



Step 1

Insert the test extraction tube into the workstation provided in the kit. Make sure that the tube is standing upright and reaches the bottom of the workstation.



Step 2

Tear off the sealing film on the extraction tube gently to avoid spilling out the liquid.



Step 3

Insert the swab into the extraction tube which contains the extraction buffer (approximately 0.3 mL).




Step 4

Roll the swab at least 6 times while pressing the head against the bottom and side of the extraction tube.



Step 5

Leave the swab in the extraction tube for 1 minute.



Step 6


Squeeze the tube several times from the outside to immerse the swab. Remove the swab.



Step 7

Insert the tip into the extraction tube tightly.

Test Procedure & Interpretation of Results



Step 1

Remove the test device from the pouch and lay it on a flat clean surface. Add 4 drops of the solution into the sample well of the test device by gently squeezing the tube.



Step 2

Start timer and wait 15 minutes. Do not read the test result after 20 minutes.


Step 3

Read your results
(See interpretation of results section).

Positive

COVID-19 Positive	The presence of two lines as control line (C) and CoV test line.
Flu A Positive	The presence of two lines as control line (C) and A test line.
Flu B Positive	The presence of two lines as control line (C) and B test line.
Flu A+B Positive	The presence of three lines as control line (C), A test line and B test line.
COVID-19 +Flu A Positive	The presence of three lines as control line (C), CoV test line and A test line.
COVID-19 +Flu B Positive	The presence of three lines as control line (C), CoV test line and B test line.
COVID-19 +Flu A & B Positive	The presence of Four lines as control line (C), CoV test line, A test and B test line.
RSV Positive	The presence of two lines as control line (C) and a (T) test line.

For COVID-19&Flu A/B




Positive

Negative

Invalid


For RSV



Positive



Negative



Invalid

Negative

The presence of only control band (C) within the result window indicates a negative result.

Invalid

If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. It is recommended that the specimen be re-tested using a new test.

Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
COVID-19/Flu A&B/RSV Ag Combo Rapid Test Cassette (Swab)	Nasal Swab	GCFCR-T525a√	Cassette	20 Tests/Kit

√ CE Marked

Zhejiang Orient Gene Biotech Co., Ltd
Address: 3787#, East Yangguang Avenue, Dipu Street,
Anji 313300, Huzhou, Zhejiang, China
Tel: +86-572-5226111
Fax: +86-572-5226222
Website: www.orientgene.com