

ToRCH (Rubella, CMV, Toxo & HSV)

The **ToRCH Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of IgG and/or IgM antibodies to Rubella virus (Rubella), Cytomegalovirus (CMV), *Toxoplasma gondii* (Toxo) and Herpes Simplex Virus (HSV) in serum, plasma or whole blood to aid in the diagnosis of ToRCH.


Materials Provided

- Rubella, CMV, Toxo & HSV Rapid Test Cassettes
- Droppers
- Buffer
- Package insert





Test Procedure

Rubella, CMV, HSV, TOXO IgM & TOXO IgG




TOXO IgM/IgG





ToRCH Combo



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Cytomegalovirus IgG Rapid Test Cassette	S/P	GCCMV(IgG)-302a	Cassette	25 Tests/Kit
Cytomegalovirus IgG Rapid Test Cassette	WB/S/P	GCCMV(IgG)-402a	Cassette	25 Tests/Kit
Cytomegalovirus IgM Rapid Test Cassette	S/P	GCCMV(IgM)-302a	Cassette	25 Tests/Kit
Cytomegalovirus IgM Rapid Test Cassette	WB/S/P	GCCMV(IgM)-402a	Cassette	25 Tests/Kit
HSV IgG Rapid Test Cassette	S/P	GCHSV(IgG)-302a	Cassette	25 Tests/Kit
HSV IgG Rapid Test Cassette	WB/S/P	GCHSV(IgG)-402a√	Cassette	25 Tests/Kit
HSV IgM Rapid Test Cassette	S/P	GCHSV(IgM)-302a√	Cassette	25 Tests/Kit
HSV IgM Rapid Test Cassette	WB/S/P	GCHSV(IgM)-402a	Cassette	25 Tests/Kit
Rubella IgG Rapid Test Cassette	S/P	GCRUB(IgG)-302a	Cassette	25 Tests/Kit
Rubella IgG Rapid Test Cassette	WB/S/P	GCRUB(IgG)-402a	Cassette	25 Tests/Kit
Rubella IgM Rapid Test Cassette	S/P	GCRUB(IgM)-302a	Cassette	25 Tests/Kit
Rubella IgM Rapid Test Cassette	WB/S/P	GCRUB(IgM)-402a	Cassette	25 Tests/Kit
Toxo IgG Rapid Test Cassette	S/P	GCTOX(IgG)-302a√	Cassette	25 Tests/Kit
Toxo IgG Rapid Test Cassette	WB/S/P	GCTOX(IgG)-402a	Cassette	25 Tests/Kit
Toxo IgM Rapid Test Cassette	S/P	GCTOX(IgM)-302a√	Cassette	25 Tests/Kit
Toxo IgM Rapid Test Cassette	WB/S/P	GCTOX(IgM)-402a	Cassette	25 Tests/Kit
Toxo IgG/IgM Rapid Test Cassette	S/P	GCTOX-302b	Cassette	25 Tests/Kit
Toxo IgG/IgM Rapid Test Cassette	S/P	GCTOX(IgG/IgM)-302a√	Cassette	20 Tests/Kit
Toxo IgG/IgM Rapid Test Cassette	WB/S/P	GCTOX-402b	Cassette	25 Tests/Kit
ToRCH IgG Combo Rapid Test Cassette	S/P	GCTOG-345a	Cassette	20 Tests/Kit
ToRCH IgM Combo Rapid Test Cassette	S/P	GCTOM-345a	Cassette	20 Tests/Kit

√ CE Marked