

Syphilis

The **Syphilis Antibody Rapid Test** is a rapid, serological, immunochromatographic assay for the qualitative detection of antibodies (IgG, IgM and IgA) to *Treponema pallidum* (TP) in human serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with TP.


Materials Provided

- Syphilis Rapid Test Strips/Cassettes
- Droppers
- Buffer
- Package insert




Test Procedure and Interpretation


Strip (S/P)



Strip (WB/S/P)




Cassette





Interpretation of Results



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Syphilis Ab Rapid Test Strip	S/P	GCSYP-301a/	Strip	50 Tests/Kit
Syphilis Ab Rapid Test Cassette	S/P	GCSYP-302a/	Cassette	25 Tests/Kit
Syphilis Ab Rapid Test Strip	WB/S/P	GCSYP-401a/	Strip	50 Tests/Kit
Syphilis Ab Rapid Test Cassette	WB/S/P	GCSYP-402a/	Cassette	25 Tests/Kit

√ CE Marked