

Gonorrhea


The **Gonorrhea Rapid Test** Cassette is intended for the qualitative detection of *Neisseria gonorrhoeae* in female cervical swab and male urethral swab specimens, as a screening test and as an aid in the diagnosis of infection with *N. Gonorrhoeae*.

Materials Provided

- Gonorrhea Rapid Test Cassettes
- Extraction tubes
- Workstation
- Swabs
- Reagent A
- Reagent B
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
Gonorrhea Rapid Test Cassette	Swab	GCGON-502b	Cassette	1.0E+7	20 Tests/Kit