

EV71

The **EV71 IgM Rapid Test** is a sandwich lateral flow chromatographic immunoassay for the qualitative detection of IgM-class antibodies to human Enterovirus 71 (EV71) in human serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with EV71.

Materials Provided

- EV71 IgM Rapid Test Cassettes
- Droppers
- Extraction tubes
- Package insert




Test Procedure and Interpretation


10 μ L serum/plasma sample
or 20 μ L whole blood sample



Stagger



2~3 drops of the diluted specimen



15 min



Positive



Negative



Invalid

Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
EV71 IgM Rapid Test Cassette	S/P	GCEV71(IgM)-302a	Cassette	25 Tests/Kit
EV71 IgM Rapid Test Cassette	WB/S/P	GCEV71(IgM)-402a	Cassette	25 Tests/Kit