

Tuberculosis


The **Tuberculosis IgG/IgM Rapid Test** is a membrane-based screening test for the rapid detection of IgM anti-Mycobacterium Tuberculosis and IgG anti-Mycobacterium Tuberculosis in human in serum, plasma or whole blood. This innovative rapid screening test is based on lateral flow immunochromatography and is among the easiest point of care (POC) assay diagnostics.

Materials Provided

- Tuberculosis Rapid Test Cassettes
- Droppers
- Buffer
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Tuberculosis IgG/IgM Rapid Test Cassette	S/P	GCTB-302a/	Cassette	25 Tests/Kit
Tuberculosis IgG/IgM Rapid Test Cassette	WB/S/P	GCTB-402a/	Cassette	25 Tests/Kit

√ CE Marked