

# CEA


The **Carcinoembryonic Antigen Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of CEA in serum, plasma or whole blood. It is intended as an aid in the assessment and diagnosis of colorectal cancer and other types of cancers, such as pancreatic, stomach, breast, lung, and certain types of thyroid and ovarian cancer.

## Materials Provided

- Carcinoembryonic Antigen Rapid Test Strips/Cassettes
- Buffer
- Droppers
- Package insert



## Test Procedure and Interpretation



## Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
Carcinoembryonic Antigen Rapid Test Strip	S/P	GECEA-301a	Strip	5 ng/mL	50 Tests/Kit
Carcinoembryonic Antigen Rapid Test Cassette	S/P	GECEA-302a	Cassette	5 ng/mL	25 Tests/Kit
Carcinoembryonic Antigen Rapid Test Strip	WB/S/P	GECEA-401a	Strip	5 ng/mL	50 Tests/Kit
Carcinoembryonic Antigen Rapid Test Cassette	WB/S/P	GECEA-402a/	Cassette	5 ng/mL	25 Tests/Kit

Zhejiang Orient Gene Biotech Co., Ltd.  
 Add: 3787#, East Yangguang Avenue, Dipu Street,  
 Anji 313300, Huzhou, Zhejiang, China  
 Tel: +86-572-5226111  
 Fax: +86-572-5226222  
 Website: www.orientgene.com

√ CE Marked