

# HEV

The **HEV IgM Rapid Test** is a lateral flow chromatographic immunoassay for the qualitative detection of IgM antibody to hepatitis E virus (HEV) in human serum or plasma. It is intended to be used as a screening test and as an aid in the diagnosis of infection with HEV.

## Materials Provided


- HEV Rapid Test Cassettes
- Droppers
- Sampling loops
- Buffer
- Package insert



## Test Procedure and Interpretation

Dip the sampling loop into specimen totally and rotate 2-3 cycles .

Put the sampling loop into the sample well (about 3μL) till touching sample pad .



Positive



Negative



Invalid

## Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
HEV Rapid Test Cassette	S/P	GCHEV-302a√	Cassette	25 Tests/Kit

√ CE Marked