

# H. pylori

The **H. pylori Rapid Test** is a sandwich lateral flow chromatographic immunoassay for the qualitative detection of antibodies (IgG, IgM, and IgA) anti-Helicobacter pylori (*H. pylori*) in human serum, plasma or whole blood and antigen in human feces. It is intended to be used as a screening test and as an aid in the diagnosis of infection with *H. pylori*.

## Materials Provided


- *H. pylori* Rapid Test Strips/Cassettes
- Droppers (For *H. pylori* Ab test use only)
- Buffer
- Package insert




## Test Procedure and Interpretation

### For *H. pylori* Ab


#### Strip



#### Cassette



## For *H. pylori* Ag



## Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
<i>H. pylori</i> Ab Rapid Test Strip	S/P	GCHP-301a√	Strip	50 Tests/Kit
<i>H. pylori</i> Ab Rapid Test Cassette	S/P	GCHP-302a√	Cassette	25 Tests/Kit
<i>H. pylori</i> Ab Rapid Test Strip	WB/S/P	GCHP-401a√	Strip	50 Tests/Kit
<i>H. pylori</i> Ab Rapid Test Cassette	WB/S/P	GCHP-402a√	Cassette	25 Tests/Kit
<i>H. pylori</i> Ag Rapid Test Strip	Feces	GCHP-601a√	Strip	25 Tests/Kit
<i>H. pylori</i> Ag Rapid Test Strip	Feces	GCHP-601Ca	Strip	25 Tests/Kit
<i>H. pylori</i> Ag Rapid Test Cassette	Feces	GCHP-602a√	Cassette	20 Tests/Kit
<i>H. pylori</i> Ag Rapid Test Cassette	Feces	GCHP-602Ca	Cassette	20 Tests/Kit

√ CE Marked