

Leishmania

The **Leishmania Ab Rapid Test** is a lateral flow immunoassay for the qualitative detection of antibodies including IgG and IgM to the subspecies of the *Leishmania donovani* (*L. donovani*), the Visceral leishmaniasis causative protozoans, in human serum or plasma. This test is intended to be used as a screening test and as an aid in the diagnosis of the disease of Visceral leishmaniasis.


Materials Provided

- Leishmania Rapid Test Strips/Cassettes
- Droppers
- Buffer
- Package insert




Test Procedure and Interpretation

Strip



Cassette



Interpretation of Results





Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Leishmania Ab Rapid Test Strip	S/P	GCKal-301a	Strip	50 Tests/Kit
Leishmania Ab Rapid Test Cassette	S/P	GCKal-302a	Cassette	25 Tests/Kit
Leishmania Ab Rapid Test Strip	WB/S/P	GCKal-401a	Strip	50 Tests/Kit
Leishmania Ab Rapid Test Cassette	WB/S/P	GCKal-402a	Cassette	25 Tests/Kit