

AFP


The **Alpha-Fetoprotein Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of AFP in serum, plasma or whole blood to aid in the diagnosis of Hepatocellular Carcinoma or fetal open neural tube defects.

Materials Provided

- Alpha-Fetoprotein Rapid Test Strips/Cassettes
- Buffer
- Droppers
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
Alpha-Fetoprotein Rapid Test Strip	S/P	GEAFP-301a	Strip	20 ng/mL	50 Tests/Kit
Alpha-Fetoprotein Rapid Test Cassette	S/P	GEAFP-302a√	Cassette	20 ng/mL	25 Tests/Kit
Alpha-Fetoprotein Rapid Test Strip	WB/S/P	GEAFP-401a	Strip	20 ng/mL	50 Tests/Kit
Alpha-Fetoprotein Rapid Test Cassette	WB/S/P	GEAFP-402a	Cassette	20 ng/mL	25 Tests/Kit

√ CE Marked

Zhejiang Orient Gene Biotech Co., Ltd.
 Add: 3787#, East Yangguang Avenue, Dipu Street,
 Anji 313300, Huzhou, Zhejiang, China
 Tel: +86-572-5226111
 Fax: +86-572-5226222
 Website: www.orientgene.com