



Micro-Albumin


The **One Step Micro-Albumin Test** is a rapid visual immunoassay for the qualitative presumptive detection of microalbumin in human urine specimens. This kit is intended to be used as an aid in the diagnosis of renal dysfunction.

Materials Provided

- One Step Micro-Albumin Test Strips
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
One Step Micro-Albumin Test Strip	Urine	GEALB-101a	Strip	20 ng/mL	100 Tests/Kit
One Step Micro-Albumin Test Strip	Urine	GEALB-101b	Strip	50 ng/mL	100 Tests/Kit