



V. cholerae


The ***V. cholerae* Antigen Rapid Test** is a rapid visual immunoassay for the qualitative, presumptive detection of *Vibrio cholerae* O1 and/or O139 antigen in human fecal specimens. This kit is intended for use as an aid in the diagnosis of *Vibrio cholerae* O1 and/or O139 infection.

Materials Provided

- *V. cholerae* Rapid Test Cassettes
- Extraction buffers
- Package insert




Test Procedure and Interpretation



V. cholerae O1



V. cholerae O1/O139





Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
<i>V. cholerae</i> O1 Ag Rapid Test Cassette	Feces	GEVCH(O1)-602a√	Cassette	20 Tests/Kit
<i>V. cholerae</i> O1/O139 Ag Combo Rapid Test Cassette	Feces	GEVCH(O1/O9)-602a√	Cassette	20 Tests/Kit

√ CE Marked