

Rapid COVID-19 Antigen Self-Test

Benefits

- Allow the individuals to collect specimens by themselves or with the assistance of others.
- Can detect all suspected case specimens or general screening specimens.
- Simple test procedure using a nasal swab from the front of the nasal cavity instead of the nasopharynx.
- With single use components, safe & easy operation and reduces the risk of pollution.
- Receive highly accurate results within 15 minutes.



Kit Contents



Swab inside sealed wrapper
(×1/2/3/5/10/15/20/25)



Test cassette in sealed pouch
(×1/2/3/5/10/15/20/25)



Plastic waste bag
(×1/2/3/5/10/15/20/25)



Extraction tube with buffer
(×1/2/3/5/10/15/20/25)




Nozzle cap
(×1/2/3/5/10/15/20/25)




Extraction tube holder
(×1/2 reusable if applicable)

Read Results




Positive result

The presence of two lines, a control line (C) and a test line (T), **even a faint line**, within the result window indicates the test is positive.



Negative result

The presence of only a control line (C) within the result window indicates the test is negative.



Void result

If the control line (C) is not visible within the result window after performing the test, the result is considered void. If you get a void result, you must perform a new test.

Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H10GE ✓	Cassette	1 Test/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H20GE ✓	Cassette	2 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H30GE ✓	Cassette	3 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H50GE ✓	Cassette	5 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H100GE ✓	Cassette	10 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H150GE ✓	Cassette	15 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H200GE ✓	Cassette	20 Tests/Kit
Rapid COVID-19 Antigen Self-Test	Nasal Swab	GCCOV-502a-H250GE ✓	Cassette	25 Tests/Kit

✓ CE Marked



Zhejiang Orient Gene Biotech Co., Ltd
Add: 3787#, East Yangguang Avenue, Dipu Street,
Anji, Huzhou, Zhejiang, China
Tel: +86-572-5303755/5303756
Fax: +86-572-5226222
Web: <http://www.orientgene.com>