

Chikungunya


The **Chikungunya IgM Rapid Test** is a lateral flow chromatographic immunoassay for the qualitative detection of IgM anti-chikungunya virus (CHIK) in human serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with CHIK.

Materials Provided

- Chikungunya IgM Rapid Test Cassettes
- Droppers
- Buffer
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Chikungunya IgM Rapid Test Cassette	S/P	GCCHK(IgM)-302a	Cassette	25 Tests/Kit
Chikungunya IgM Rapid Test Cassette	WB/S/P	GCCHK(IgM)-402a√	Cassette	25 Tests/Kit

√ CE Marked

Zhejiang Orient Gene Biotech Co., Ltd.
 Add: 3787#, East Yangguang Avenue, Dipu Street,
 Anji 313300, Huzhou, Zhejiang, China
 Tel: +86-572-5226111
 Fax: +86-572-5226222
 Website: www.orientgene.com