

Strep A

The **Strep A Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of Strep A antigen from throat swab specimens to aid in the diagnosis of Group A Streptococcal infection.


Materials Provided

- Strep A Rapid Test Strips/Cassettes
- Sterile swabs
- Extraction tubes and dropper tips
- Reagent A
- Reagent B
- Negative control (when applicable)
- Positive control (when applicable)
- Workstation
- Package insert



Test Procedure and Interpretation

Test Procedure



Interpretation of Results





Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Strep A Rapid Test Strip	Throat Swab	GCSTR-501a√	Strip	25 Tests/Kit
Strep A Rapid Test Strip	Throat Swab	GCSTR-501Ca√	Strip	25 Tests/Kit
Strep A Rapid Test Cassette	Throat Swab	GCSTR-502a√	Cassette	20 Tests/Kit
Strep A Rapid Test Cassette	Throat Swab	GCSTR-502Ca√	Cassette	20 Tests/Kit

√ CE Marked