

Chagas


The **Chagas Antibody Rapid Test** is a lateral flow chromatographic immunoassay for the qualitative detection of IgG anti-*Trypanosoma cruzi* (*T. cruzi*) in human serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with *T. Cruzii*.

Materials Provided

- Chagas Antibody Rapid Test Cassettes
- Droppers
- Buffer
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Chagas Ab Rapid Cassette	S/P	GCCHA-302a	Cassette	25 Tests/Kit
Chagas Ab Rapid Cassette	WB/S/P	GCCHA-402a√	Cassette	25 Tests/Kit

√ CE Marked