

RF


The **Rheumatoid Factor IgM Rapid Test** is a lateral flow chromatographic immunoassay for the qualitative detection of IgM antibodies in human serum to Rheumatoid Factor (RF) antigen. It is intended to be used as a screening test and as an aid in the diagnosis of rheumatoid arthritis.

Materials Provided

- Rheumatoid Factor IgM Rapid Test Cassettes
- Buffer
- Droppers
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Rheumatoid Factor IgM Rapid Test Cassette	S/P	GGRF(IgM)-302a	Cassette	25 Tests/Kit