

HBsAg/HCV/HIV/Syphilis Combo


The **HBsAg/HCV/HIV/Syphilis Combo Rapid Test** is a rapid, immunochromatographic assay for the qualitative detection of HBsAg, HCV antibodies, HIV1/2 virus-antibodies and/or *Treponema pallidum* antibodies respectively in human whole blood, serum or plasma samples. It is intended for professional use as an aid in diagnosing HBV, HCV, HIV and/or Syphilis infections.

Materials Provided

- HBsAg/HCV/HIV/Syphilis Combo Rapid Test Cassettes
- Droppers
- Buffer
- Package insert



Test Procedure and Interpretation



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
HBsAg/HCV Combo Rapid Test Cassette	WB/S/P	GCHBC-402a	Cassette	25 Tests/Kit
HCV/HIV Combo Rapid Test Cassette	WB/S/P	GCHCI-402a	Cassette	25 Tests/Kit
HIV/HCV/HBsAg/Syphilis Combo Rapid Test Cassette	S/P	GCHBCISY-345a	Cassette	20 Tests/Kit
HIV/HCV/HBsAg/Syphilis Combo Rapid Test Cassette	WB/S/P	GCHBCISY-445a	Cassette	20 Tests/Kit