

HBV

The **Hepatitis B Virus Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of Hepatitis B Surface Antigen (HBsAg), Hepatitis B Surface Antibody (HBsAb), Hepatitis B Envelope Antigen (HBeAg), Hepatitis B Envelope Antibody (HBeAb) and Hepatitis B core Antibody (HBcAb) in serum, plasma or whole blood.


Materials Provided

- HBcAb/HBeAb/HBeAg/HBsAb/ HBsAg Rapid Test Strips/Cassettes
- Droppers
- Buffer (For HBsAb Cassette WB/S/P, HBsAg Strip /Cassette WB/S/P & Combo test use only)
- Package insert




Test Procedure and Interpretation


HBcAb/HBeAb/HBeAg Test




HBsAb strip




HBsAb cassette




HBsAg strip (S/P)




HBsAg strip (WB/S/P)




HBsAg cassette



HBV Combo



Interpretation of Results:



Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
One Step HBcAb Test Strip	S/P	GCHBcb-302a	Cassette	2 NCU	25 Tests/Kit
One Step HBcAb Test Cassette	S/P	GCHBcb-302b	Cassette	8 NCU	25 Tests/Kit
One Step HBeAb Test Cassette	S/P	GCHBeb-302a	Cassette	2 NCU	25 Tests/Kit
One Step HBeAb Test Cassette	S/P	GCHBeb-302b	Cassette	8 NCU	25 Tests/Kit
One Step HBeAg Test Cassette	S/P	GCHBeg-302a	Cassette	0.5 NCU	25 Tests/Kit
One Step HBsAb Test Strip	S/P	GCHBsb-301a	Strip	30 mIU/mL	50 Tests/Kit
One Step HBsAb Test Cassette	S/P	GCHBsb-302a	Cassette	30 mIU/mL	25 Tests/Kit
HBsAb Rapid Test Strip	WB/S/P	GCHBsb-401a	Strip	30 mIU/mL	50 Tests/Kit
HBsAb Rapid Test Cassette	WB/S/P	GCHBsb-402a	Cassette	30 mIU/mL	25 Tests/Kit
HBsAb Rapid Test Cassette	WB/S/P	GCHBsb-402b	Cassette	20 mIU/mL	25 Tests/Kit
One Step HBsAg Test Strip	S/P	GCHBsg-301a	Strip	1 ng/mL	50 Tests/Kit
One Step HBsAg Test Cassette	S/P	GCHBsg-302a	Cassette	1 ng/mL	25 Tests/Kit
HBsAg Rapid Test Strip	WB/S/P	GCHBsg-401a	Strip	1 ng/mL	50 Tests/Kit
HBsAg Rapid Test Cassette	WB/S/P	GCHBsg-402a	Cassette	1 ng/mL	25 Tests/Kit