

PSA

The **Prostate Specific Antigen Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of Prostate Specific Antigen (PSA) in human serum, plasma or whole blood. It is intended for use at point of care facilities to measure total prostate specific antigen (tPSA).

Materials Provided


- Prostate Specific Antigen Rapid Test Strips/Cassettes
- Buffer
- Droppers
- Package insert




Test Procedure and Interpretation

Qualitative Test


S/P Strip Test




WB/S/P Strip Test



Cassette Test



Semi-Quantitative Test



Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
Prostate Specific Antigen Rapid Test Strip	S/P	GEPSA-301a√	Strip	4 ng/mL	50 Tests/Kit
Prostate Specific Antigen Rapid Test Cassette	S/P	GEPSA-302a√	Cassette	4 ng/mL	25 Tests/Kit
Prostate Specific Antigen Rapid Test Strip	WB/S/P	GEPSA-401a√	Strip	4 ng/mL	50 Tests/Kit
Prostate Specific Antigen Rapid Test Cassette	WB/S/P	GEPSA-402a√	Cassette	4 ng/mL	25 Tests/Kit
Prostate Specific Antigen Semi-Quantitative Rapid Test Cassette	S/P	GEPSA-302b	Cassette	4 ng/mL, 10 ng/mL	25 Tests/Kit
Prostate Specific Antigen Semi-Quantitative Rapid Test Cassette	WB/S/P	GEPSA-402b	Cassette	4 ng/mL, 10 ng/mL	25 Tests/Kit

√ CE Marked