

Dengue

The **Dengue Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of IgG, IgM antibodies and/or NS1 antigen to dengue virus in human whole blood, serum or plasma as an aid in the diagnosis of dengue infections.


Materials Provided

- Dengue Rapid Test Cassettes
- Droppers
- Buffer
- Package insert




Test Procedure and Interpretation


Dengue NS1




Dengue IgM/IgG




Dengue NS1+IgG/IgM






Interpretation of Results

IgG/IgM Antibody Test:



NS1 Antigen Test:



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
Dengue IgM/IgG Rapid Test Cassette	WB/S/P	GC DEN(ab)-402c√	Cassette	25 Tests/Kit
Dengue NS1 Ag Rapid Test Cassette	WB/S/P	GC DEN(NS)-402c√	Cassette	25 Tests/Kit
Dengue NS1 and IgM/IgG Combo Rapid Test Cassette	WB/S/P	GC DEN-425a√	Cassette	20 Tests/Kit

√ CE Marked