

Utilize a series of visual, semi-quantitative tests simultaneously to analyze urine properties that aid in detecting/monitoring conditions such as:

- Diabetes
- Metabolic Abnormalities
- Liver Disease
- Kidney Function
- Urinary Tract Infection

Results read within 2 minutes. These tests can be used along with other diagnostic information to rule out certain disease states and determine if further microscropic analysis is needed.

## **Testing Parameters**

Urobilinogen	Protein
Bilirubin	Nitrite
Ketone	Ascorbi
Blood	Leukocy

Glucose Specific Gravity ic Acid pH

ytes

## Urinalysis Tests are now available Conveniently Dropship from Houston



510(k) cleared, CLIA Waived Prescription Use Only

