

# FOB & TF


The **Fecal Occult Blood/Transferrin Rapid Test** is a one step chromatographic immunoassay for the qualitative determination of human hemoglobin (Hb) or human transferrin (Tf) in stool samples to detect gastrointestinal bleeding. The test is intended for use by health care professionals as an aid in the diagnosis of premature syndrome linked to the intestinal bleeding.

## Materials Provided


- FOB/Transferrin Rapid Test Strips/Cassettes
- Buffers
- Droppers
- Package insert




## Test Procedure and Interpretation



### FOB or TF Test



### FOB and TF Combo





## Ordering Information

| Product Description                                          | Specimen | Catalog No.    | Format   | Cut-off Value | Kit Size     |
|--------------------------------------------------------------|----------|----------------|----------|---------------|--------------|
| Fecal Occult Blood Rapid Test Strip                          | Feces    | GEFOB-601b√†   | Strip    | 50 ng/mL      | 25 Tests/Kit |
| Fecal Occult Blood Rapid Test Strip                          | Feces    | GEFOB-601c√    | Strip    | 100 ng/mL     | 25 Tests/Kit |
| Fecal Occult Blood Rapid Test Strip                          | Feces    | GEFOB-601d     | Strip    | 200 ng/mL     | 25 Tests/Kit |
| Fecal Occult Blood Rapid Test Cassette                       | Feces    | GEFOB-602b√†   | Cassette | 50 ng/mL      | 20 Tests/Kit |
| Fecal Occult Blood Rapid Test Cassette                       | Feces    | GEFOB-602c√    | Cassette | 100 ng/mL     | 20 Tests/Kit |
| Fecal Occult Blood Rapid Test Cassette                       | Feces    | GEFOB-602d     | Cassette | 200 ng/mL     | 20 Tests/Kit |
| Fecal Occult Blood Rapid Test Cassette                       | Feces    | GEFOB-602h     | Cassette | 150 ng/mL     | 20 Tests/Kit |
| Transferrin Rapid Test Strip                                 | Feces    | GETF-601a√     | Strip    | 10 ng/mL      | 25 Tests/Kit |
| Transferrin Rapid Test Cassette                              | Feces    | GETF-602a√     | Cassette | 10 ng/mL      | 20 Tests/Kit |
| Fecal Occult Blood and Transferrin Combo Rapid Test Cassette | Feces    | GEFOB/TF-602a√ | Cassette | 50/10 ng/mL   | 20 Tests/Kit |

√ CE Marked    † Cleared for US 510(k)