

# Chlamydia

The **Chlamydia Trachomatis Antigen Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of Chlamydia in female cervical swab, male urethral swab and male urine specimens. The product can detect the Chlamydia serovars (D, E, F, H, I, K, G, J) and intended as a screening test and as an aid in the diagnosis of Chlamydia infection.


## Materials Provided

- Chlamydia Ag Rapid Test Cassettes
- Extraction tubes
- Swabs
- Workstation
- Extraction buffer A
- Extraction buffer B
- Package insert



## Test Procedure and Interpretation

For swab specimen



For male urine specimen, please refer to the package insert carefully.

## Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
Chlamydia Ag Rapid Test Cassette	Swab/Urine	GCCHL-502a/	Cassette	$4.8 \times 10^3$ IFU/mL	20 Tests/Kit
Zhejiang Orient Gene Biotech Co., Ltd. Add: 3787#, East Yangguang Avenue, Dipu Street, Anji 313300, Huzhou, Zhejiang, China Tel: +86-572-5226111 Fax: +86-572-5226222 Website: www.orientgene.com					√ CE Marked