

HIV

The **Human Immunodeficiency Virus Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of HIV p24 antigen and/or antibodies (IgG, IgM and IgA) to HIV-1 and/or HIV-2 and/or Subtype O in serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with HIV.


Materials Provided

- HIV Rapid Test Strips/Cassettes
- Droppers
- Buffer
- Package insert




Test Procedure and Interpretation


Strip





Cassette





Interpretation of Results



Tri-line



HIV Ag/Ab



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
HIV 1/2 Ab Rapid Test Strip	S/P	GCHIV-301a	Strip	50 Tests/Kit
HIV 1/2 Ab Rapid Test Cassette	S/P	GCHIV-302a/	Cassette	25 Tests/Kit
HIV 1/2 Ab Rapid Test Strip	WB/S/P	GCHIV-401a	Strip	50 Tests/Kit
HIV 1/2 Ab Rapid Test Cassette	WB/S/P	GCHIV-402a/	Cassette	25 Tests/Kit
HIV 1/2 Ab Tri-line Rapid Test Cassette	WB/S/P	GCHIV-GT402a	Cassette	25 Tests/Kit
HIV 1/2/O Ab Tri-line Rapid Test Cassette	S/P	GCHIV-T302b	Cassette	25 Tests/Kit
HIV 1/2/O Ab Tri-line Rapid Test Cassette	WB/S/P	GCHIV-T402a	Cassette	25 Tests/Kit
HIV Ab/Ag Rapid Test Cassette	WB/S/P	GCHIV(Ag/Ab)-402a	Cassette	25 Tests/Kit

√ CE Marked