

HCV

The **HCV Antibody Rapid Test** is a sandwich lateral flow chromatographic immunoassay for the qualitative detection of antibodies (IgG, IgM, and IgA) anti-Hepatitis C Virus (HCV) in human serum, plasma or whole blood. It is intended to be used as a screening test and as an aid in the diagnosis of infection with HCV.


Materials Provided

- HCV Rapid Test Strips/Cassettes
- Droppers
- Buffer
- Package insert




Test Procedure and Interpretation

Strip



Cassette



Interpretation of Results



Ordering Information

Product Description	Specimen	Catalog No.	Format	Kit Size
HCV Ab Rapid Test Strip	S/P	GCHCV-301a	Strip	50 Tests/Kit
HCV Ab Rapid Test Cassette	S/P	GCHCV-302a/	Cassette	25 Tests/Kit
HCV Ab Rapid Test Strip	WB/S/P	GCHCV-401a	Strip	50 Tests/Kit
HCV Ab Rapid Test Cassette	WB/S/P	GCHCV-402a/	Cassette	25 Tests/Kit

√ CE Marked

Zhejiang Orient Gene Biotech Co., Ltd.
 Add: 3787#, East Yangguang Avenue, Dipu Street,
 Anji 313300, Huzhou, Zhejiang, China
 Tel: +86-572-5226111
 Fax: +86-572-5226222
 Website: www.orientgene.com