



Hepatitis B Virus Detection Kit (Fluorescence PCR)

Intended use

This kit can be used for the quantitative detection of Hepatitis B virus DNA (genotypes A-H) in human plasma or serum. The results are only for clinical reference and should not be used as conclusive evidence for the diagnosis or exclusion of diseases.

Features:



- High Sensitivity
- High Specificity
- Reliable Result

Specification

- | | |
|---|--|
| • Protocol Duration:
90-120 minutes | • Analytical Sensitivity (LoD):
20 IU/mL |
| • Storage:
-25°C to -15°C | • Shelf Life:
12 months from the date of manufacture |



Product Description	Specimen	Catalog No.	Format	Kit Size
Hepatitis B Virus Detection Kit (Fluorescence PCR)	Plasma (EDTA anticoagulant) or serum	101A0226EY	Kit	48 tests/kit



Wide range of applicable models:

It is suitable for ABI, Roche, Bio-Rad, Bioer, Hongshi, Molarray and similar multi-channel fluorescent PCR.

Equipment on sale

HG-P960 Real-time PCR system

- Automatic pop-up sample bin
- Intelligent adjustable hot cover
- 6 partition thermal cycling module
- Full adaptable software system
- Top imaging photoelectric detection



OG-P100 Real-time PCR system

- General consumables matching, easy to use
- No need to calibrate regularly
- Fast mode can complete the test in 20 minutes
- 10 inch screen is easy to operate and save space
- Self owned special chip to optimize instrument structure
- Independent research and development, flexible combination and customization

HG-P320 Real-time PCR system

- The experimental results can be exported directly.
- 4.7-inch high-definition TFT color touch screen, and embedded operating system.
- 4 channels and double 16-well blocks design, can run two different programs at the same time.
- Powerful software analysis function, which can be used for Quantitative Analysis, Melting Curve Analysis, etc.

