

## *M. pneumoniae*


The **Mycoplasma IgM Rapid Test** is a rapid chromatographic immunoassay for the qualitative detection of IgM antibodies to *Mycoplasma pneumoniae* antigen in serum or plasma to aid in the diagnosis of recent *Mycoplasma pneumoniae* infection.

### Materials Provided

- Mycoplasma IgM Rapid Test Cassettes
- Droppers
- Buffer
- Package insert



### Test Procedure and Interpretation



### Ordering Information

| Product Description                          | Specimen | Catalog No.     | Format   | Kit Size     |
|----------------------------------------------|----------|-----------------|----------|--------------|
| <i>M. pneumoniae</i> IgM Rapid Test Cassette | S/P      | GCMP(IgM)-302a/ | Cassette | 25 Tests/Kit |
|                                              |          |                 |          | √ CE Marked  |