

Influenza

The **Influenza Antigen Rapid Test** is an *in vitro* immunochromatographic assay for the qualitative detection of influenza A (including the subtype H1N1) and/or B nucleoprotein antigens in nasopharyngeal (NP) swab, nasal swab, and nasal wash/aspirate specimens. It is intended to aid in the rapid differential diagnosis of influenza A and/or B viral infections.

Materials Provided

- Influenza A/A&B Rapid Test Strips/Cassettes
- Sterile swabs
- Extraction tubes
- Dropper tips (for cassette use only)
- Positive control swab (when applicable)
- Negative control swab (when applicable)
- Workstation
- Buffer
- Package insert



Test Procedure and Interpretation

Influenza A



Influenza A&B





Ordering Information

Product Description	Specimen	Catalog No.	Format	Cut-off Value	Kit Size
Influenza A Ag Rapid Test Strip	Nasal/Throat Swabs	GCFLU(A)-501a√	Strip	1.5×10^4 TCID ₅₀	25 Tests/Kit
Influenza A Ag Rapid Test Cassette	Nasal/Throat Swabs	GCFLU(A)-502a√	Cassette	1.5×10^4 TCID ₅₀	20 Tests/Kit
Influenza A&B Ag Rapid Test Strip	Nasal/Throat Swabs	GCFLU(A/B)-501a√	Strip	1.5×10^4 TCID ₅₀ / 1.5×10^5 TCID ₅₀	25 Tests/Kit
Influenza A&B Ag Rapid Test Cassette	Nasal/Throat Swabs	GCFLU(A/B)-502a√	Cassette	1.5×10^4 TCID ₅₀ / 1.5×10^5 TCID ₅₀	20 Tests/Kit
Influenza A&B Ag Rapid Test Cassette	Nasal/Throat Swabs	GCFLU(A/B)-502Ca√	Cassette	1.5×10^4 TCID ₅₀ / 1.5×10^5 TCID ₅₀	20 Tests/Kit

√ CE Marked