



Urinalysis Strips



Testing Parameters: LEU | NIT | URO | PRO | PH | BLD
SG | KET | BIL | GLU | VC | MA | CRE | CA
Featured Combinations: UTI | ACR | Vaginal pH

FDA 510K cleared, CLIA Waived

Urine Analyzer



Healgen 500

- 120 T/hour or 60 T/hour
- High-luminance white LED, long life and high stability
- Built-in high-speed thermal printer
- Support external barcode scanner



Healgen 501

- Portable handheld
- Test speed: 60 T/hour
- High-luminance white LED, long life and high stability
- Built-in rechargeable lithium battery



Healgen 800

- Up to 500 tests/hour
- Large LCD color touchscreen
- Strip lockout to prevent using strips from another brand
- RS232 port for data and LIS