



A contribution to human health.




Cardiac Marker Test

Zhejiang Orient Gene Biotech Co., Ltd
Zhejiang, China
+86-572-5303755/5303756
www.orientgene.com


Product Name	Specimen	Catalog No.	Format	Package	Sensitivity	Certificate
Troponin I Test	Serum/Plasma	GDTRO-301	Strip	50T/Box	0.5ng/ml	CE
		GDTRO-302	Device	25T/Box	0.5ng/ml	CE
	Whole Blood/Serum/Plasma	GDTRO-401	Strip	50T/Box	0.5ng/ml	CE
		GDTRO-402	Device	25T/Box	0.5ng/ml	CE
Myoglobin test	Whole Blood/Serum/Plasma	GDMYO-402	Device	25T/Box	50ng/ml	CE
CK-MB test	Serum/Plasma	GDCKM-302	Device	25T/Box	5ng/ml	
	Whole Blood/Serum/Plasma	GDCKM-402	Device	25T/Box	5 ng/ml	CE
CRP Semi-Quantitative C-Reactive Protein Test	Whole Blood/Serum/Plasma	GDCRP-402	Device	25T/Box	3mg/L 10mg/L	CE
Cardiac Myoglobin /CK- MB/cTnI Combo test	Whole Blood/Serum/Plasma	GDCAR-435	Device	25T/Box	50/5/0.5ng/ml	CE

TEST PROCEDURE


Troponin I Test




CRP Test



Myoglobin



CK-MB Test



Cardiac Marker Test

